

**Precision balance Type PBS/PBJ****Kern & Sohn****Precision balance Type PBS/PBJ**

- Only PBJ: Automatic internal adjustment in case of a change in temperature or time-controlled at defined intervals
- Metal housing: robust and sturdy
- Large backlit, LCD display, digit height 14 mm
- Dosage aid: High-stability mode and other settings can be selected
- Weighing with tolerance range (checkweighing): Input of an upper/lower limit value. A visual signal assists with portion division, dispensing or grading
- Totalizing of individual weighing results
- Identification number: 4 numeric positions, printed on calibration protocol
- Automatic data output to the PC/printer every time the balance is steady
- Draft shield standard, on all models with readout $d = 0,001 \text{ g}$, weighing space (W x D x H): 180 x 193 x 87mm

Single-cell advanced technology:

- Fully automatic manufactured weighing cell from one piece of material
- Stable temperature behaviour
- Short stabilisation time: 3 sec
- Shock proof construction
- High corner load performance

Specifications**Overall dimensions, without draft shield (W x D x H): 209 x 322 x 78mm****Dimensions weighing space (W x D x H): 180 x 193 x 87mm****Permissible ambient temperature: 10°C / 30°C****Weight: approx. 3kg**

Code	Description	Packaging
LLG06245038	Precision balance Type PBS/PBJ, Type PBJ 4200-2M, 4200 g	1 pz.
LLG06255611	Precision balance Type PBS/PBJ, Type PBS 4200-2M, 4200 g	1 pz.
LLG07910128	Precision balance Type PBS/PBJ, Type PBJ 620-3M, 620 g	1 pz.
LLG06236719	Precision balance Type PBS/PBJ, Type PBJ 6200-2M, 6200 g	1 pz.
LLG09901445	Precision balance Type PBS/PBJ, Type PBS 6200-2M, 6200 g	1 pz.
LLG07657889	Precision balance Type PBS/PBJ, Type PBS 620-3M, 620 g	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

