

**Vortex-Genie-series**

Scientific Industries

Vortex-Genie®-series

Robust and reliable shakers for speedy mixing and shaking of samples. Multiple accessories: e.g. adapter tops for test tubes, beakers, centrifuge tubes etc. Adapter tops can be changed easily with the snap fit system.

Vortex-Genie® 2: Variable speed for precise mixing from gentle to vigorous. Broad range of attachments for most tubes, plates and other containers. Hands-Free or Touch On control.

Vortex-Genie® 2T with an integrated timer: Timer function: (Touch On) 1-60 seconds, (Hands-Free) 1-60 minutes or continuous. Unit automatically stops at the end of user programmed time.

Digital Vortex-Genie® 2 with digital control and display of speed and time: Digital display provides accuracy, reproducibility and repeatability. Timer function: (Touch On) 1-99 seconds, (Hands-Free) 1-99 minutes or continuous. Specially configured unit for Applied Biosystems protocol available.

Vortex Genie® Pulse: Fully programmable vortex mixer with ON/OFF feature. Pulse 0.1 sec to 99 min.



Scope of supply: Vortexer, 3-inch Platform with rubber cover, Pop-off cup

Specifications Vortex-Genie® 2 / Vortex-Genie® 2T / Digital Vortex-Genie® 2 / Vortex Genie® Pulse

Shaking movement: orbital

Speed display/Timer: Scale/No / Scale/Scale / Digital/Digital

Speed range: 600-2700rpm / 600-2700rpm / 500-2850rpm / 500-2850rpm

Timer

Interval operation:

Hands-free mode: - / 1 to 60s / 1 to 99s / 0.1 s to 99 min

- / 1 to 60min / 1 to 99min / 1 to 99 min

Dimensions (L x W x H): 165 x 122 x 165mm

Permissible ambient conditions: 0 to 38°C at 95% rh



Weight: 4kg
Power supply: 230V, 50/60Hz

Code	Description	Packaging
LLG06235684	Vortex Genie® 2-Mixer 240 V, 50 Hz with swiss plug	1 pz.
LLG06242793	Digital VORTEX GENIE 2-Mixer 230 V, 50-60 Hz, 500-2850 U/min	1 pz.
LLG06251939	Vortex Genie® 2-Mixer 230 V, 50 Hz, with UK-Plug	1 pz.
LLG06252823	Vortex Genie 2T 230 V, incl.standard attachment	1 pz.



Code	Description	Packaging
LLG06265224	Vortex Genie® 2-Mixer 240 V, 50 Hz with Australian plug	
LLG06265225	Vortex Genie® 2T 230 V, with standard attachment, with Australian plug	
LLG06265226	Digital Vortex Genie® 2-Mixer 230 V, Australian plug	
LLG06287671	Vortex Genie® Pulse 230 V, EU plug	1 pz.
LLG06287674	Vortex Genie® Pulse 230 V, UK plug	1 pz.
LLG06287675	Vortex Genie® Pulse 230 V, CH plug	1 pz.
LLG09730060	Vortex Genie 2-Mixer 230 V, 50 Hz	1 pz.
LLG09730120	Digital Vortex Genie® 2-Mixer 230 V, 50-60 Hz, 500-2850 U/min with UK plug	1 pz.
LLG06240599	Digital Vortex-Genie 2 230 V, Swiss plug	1 pz.
LLG06252946	Vortex-Genie 2T 230 V, incl. standard attachment, with UK plug	1 pz.
LLG09730121	Vortex Genie® 2T 230 V, with standard attachment, with CH plug	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

