



LLG-Buckets, PP

LLG-Buckets, PP

White. With lid and temper-evident closure.

For transport, packaging and storage of foodstuffs, building materials and non-corrosive chemicals.

Stability guaranteed when filling up with material up to 80 °C.



Code	Description	Packaging
LLG09040030	Buckets, PP, 1 Litres, Diam. 131 mm, Height 127 mm	1 pz
LLG09040032	Buckets, PP, 3 Litres, Diam. 198 mm, Height 141 mm	1 pz
LLG09040034	Buckets, PP, 10 Litres, Diam. 266 mm, Height 265 mm	1 pz
LLG09040031	Buckets, PP, 2 Litres, Diam. 179 mm, Height 133 mm	1 pz
LLG09040033	Buckets, PP, 5 Litres, Diam. 224 mm, Height 193 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

