



Beakers, low form, heavy duty, Pyrex

Beakers, low form, heavy duty, Pyrex®

- Mechanically stronger than standard beakers
- Designed for longer life and increased user safety
- Reinforced rims and spout, less prone to breakage during washing etc.
- Thicker walls and heavier bases
- Manufactured from Pyrex® borosilicate glass
- White graduations and marking spot



Code	Description	Packaging
LLG06252670	Beakers, low form, heavy duty, Pyrex®, Nominal capacity 1000 ml	6 pz.
LLG06230404	Beakers, low form, heavy duty, Pyrex®, Nominal capacity 250 ml	12 pz.
LLG06230405	Beakers, low form, heavy duty, Pyrex®, Nominal capacity 400 ml	12 pz.
LLG06235736	Beakers, low form, heavy duty, Pyrex®, Nominal capacity 600 ml	6 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

