



Flasks, Büchner, plain side-arm, Pyrex

Flasks, Büchner, plain side-arm, Pyrex®

- Complies with ISO 6556
- Büchner type flasks for vacuum filtration
- Heavy walls provide increased mechanical strength
- "Studded" base to minimise surface abrasion
- With integral glass side-arm to accept 9mm bore flexible tubing



Code	Description	Packaging
LLG09014036	Flasks, Büchner, plain side-arm, Pyrex®, 500 ml, Neck diam. 35 mm, External flask diam. 105 mm, Height 185 mm	10 pz.
LLG09014037	Flasks, Büchner, plain side-arm, Pyrex®, 1000 ml, Neck diam. 45 mm, External flask diam. 135 mm, Height 230 mm	5 pz.
LLG09014038	Flasks, Büchner, plain side-arm, Pyrex®, 3000 ml, Neck diam. 68 mm, External flask diam. 187 mm, Height 335 mm	5 pz.
LLG09014035	Flasks, Büchner, plain side-arm, Pyrex®, 250 ml, Neck diam. 35 mm, External flask diam. 85 mm, Height 155 mm	10 pz.
LLG09014039	Flasks, Büchner, plain side-arm, Pyrex®, 5000 ml, Neck diam. 68 mm, External flask diam. 220 mm, Height 385 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

