

**Digital thermo-hygrometer for room and outdoor measurement**

TFA Dostmann GmbH &amp; Co.KG

**Digital thermo-hygrometer for room and outdoor measurement**

For simultaneous monitoring of indoor and outdoor temperature/humidity.

- Measurement of outdoor data via cable sensor (cable length approx. 1.5 m)
- Display of minimum and maximum values and observation over a selectable period of time
- Min-/Max function
- For hanging and standing

**Scope of supply:** Thermo-hygrometer, battery

**Specifications**

**Temperature range:** -10 ... 60 °C

**Temperature accuracy:**  $\pm 0.8$  °C at 10 ... 30 °C, otherwise  $\pm 1$  °C

**Temperature resolution:** 0.1 °C

**Humidity range:** 10 ... 99 %

**Humidity accuracy:**  $\pm 3.5$  % at 30 ... 80 %, otherwise  $\pm 5$  %

**Humidity resolution:** 1 %

**Power supply:** 1 x 1.5 V AAA

**Dimensions (W x D x H):** 102 x 21 x 110 mm

**Weight:** 144 g



Type Digital thermometer/hygrometer



Code	Description	Packaging
LLG09726418	Digital thermo-hygrometer for room and outdoor measurement	1 pc



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

