

**Digital thermo-hygrometer for room and outdoor measurement**

TFA Dostmann GmbH & Co.KG

Digital thermo-hygrometer for room and outdoor measurement

For simultaneous monitoring of indoor and outdoor temperature/humidity.

- Measurement of outdoor data via cable sensor (cable length approx. 1.5 m)
- Display of minimum and maximum values and observation over a expeely selectable period of time
- Min-/Max function
- For hanging and standing

Scope of supply: Thermo-hygrometer, battery

Specifications

Temperature range: -10 ... 60 °C

Temperature accuracy: ± 0.8 °C at 10 ... 30 °C, otherwise ± 1 °C

Temperature resolution: 0.1 °C

Humidity range: 10 ... 99 %

Humidity accuracy: ± 3.5 % at 30 ... 80 %, otherwise ± 5 %

Humidity resolution: 1 %

Power supply: 1 x 1.5 V AAA

Dimensions (W x D x H): 102 x 21 x 110 mm

Weight: 144 g

Type Digital thermometer/hygrometer



Code	Description	Packaging
LLG09726418	Digital thermo-hygrometer for room and outdoor measurement	1 pc



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

