

**LLG-2in1 KITs with Screw Neck Vials ND8 (small opening)****LLG Labware****LLG-2in1 KITs with Screw Neck Vials ND8 (small opening)**

2in1 KITs contain 100 vials (1.5 ml) and 100 closures in an orange PP-Box. As both components are always required at the same time, 2in1 KITs are a convenient way to obtain and store all required items for analysis. All advantages of the individual components (cleanroom packaging of the vials, recloseability of the packaging, etc.) remain unchanged.

Any other combination of vial and seal as a 2in1 KIT available.



Code	Description	Packaging
LLG06223506	2in1 KITs with Screw Neck Vials ND8 (small opening), Type 3in1 Kit, clear, 1.5 ml, Cover type Screw Cap, black, hole, Silicon white /PTFE blue, slitted (Septa unassembled)	100 pz.
LLG06238965	2in1 KITs with Screw Neck Vials ND8 (small opening), Type 2in1 Kit, clear, 1.5 ml, Cover type Screw Cap, black, hole, Silicon white / PTFE red	100 pz.
LLG09003557	2in1 KITs with Screw Neck Vials ND8 (small opening), Type 2in1 Kit, amber, 1.5 ml, Cover type Screw Cap, black, hole, Silicon white / PTFE red	100 pz.
LLG09003558	2in1 KITs with Screw Neck Vials ND8 (small opening), Type 2in1 Kit, amber, labelling field, 1.5 ml, Cover type Screw Cap, black, hole, Silicon white / PTFE red	100 pz.
LLG09003559	2in1 KITs with Screw Neck Vials ND8 (small opening), Type 2in1 Kit, clear, labelling field, 1.5 ml, Cover type Screw Cap, black, hole, Silicon white / PTFE red	100 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

