

**Stirrer bearings with GL male thread, DURAN tubing****Lenz****Stirrer bearings with GL male thread, DURAN® tubing**

Made of DURAN® tubing. Other materials used: PTFE (bearing), silicone (gasket) and stirrer guides, GL male thread. Cone size 29/32. Guide bore 10mm i.d.



Code	Description	Packaging
LLG09197355	Stirrer bearings with GL male thread, DURAN® tubing, Ø agitator shaft 10 mm, Length 70 mm, Neck thread 25 GL, Cone 29/32 NS	1 pz.
LLG09197356	Stirrer bearings with GL male thread, DURAN® tubing, Ø agitator shaft 16 mm, Length 85 ± 3 mm, Neck thread 32 GL, Cone 45/40 NS	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

