




Magnetic stirring bars, PTFE, cylindrical

Magnetic stirring bars, PTFE, cylindrical

PTFE-coated. Cylindrical. Round, smooth surface. Versatile. Autoclavable.
Further sizes available on request.



Code	Description	Packaging
LLG09197500	Magnetic stirring bars, PTFE, cylindrical, Diam. 1.7 mm, Length 8 mm	
LLG09197501	Magnetic stirring bars, PTFE, cylindrical, Diam. 1.7 mm, Length 15 mm	
LLG09197507	Magnetic stirring bars, PTFE, cylindrical, Diam. 3 mm, Length 10 mm	
LLG09197505	Magnetic stirring bars, PTFE, cylindrical, Diam. 3.5 mm, Length 3.5 mm	
LLG09197512	Magnetic stirring bars, PTFE, cylindrical, Diam. 4.5 mm, Length 12 mm	
LLG09197503	Magnetic stirring bars, PTFE, cylindrical, Diam. 2 mm, Length 5 mm	
LLG09197504	Magnetic stirring bars, PTFE, cylindrical, Diam. 2 mm, Length 7 mm	
LLG09197502	Magnetic stirring bars, PTFE, cylindrical, Diam. 2.5 mm, Length 2.5 mm	
LLG09197506	Magnetic stirring bars, PTFE, cylindrical, Diam. 3 mm, Length 6 mm	
LLG09197508	Magnetic stirring bars, PTFE, cylindrical, Diam. 3 mm, Length 13 mm	
LLG09197509	Magnetic stirring bars, PTFE, cylindrical, Diam. 3 mm, Length 8 mm	
LLG09197514	Magnetic stirring bars, PTFE, cylindrical, Diam. 4.5 mm, Length 15 mm	
LLG09197510	Magnetic stirring bars, PTFE, cylindrical, Diam. 6 mm, Length 10 mm	
LLG09197515	Magnetic stirring bars, PTFE, cylindrical, Diam. 6 mm, Length 15 mm	
LLG09197520	Magnetic stirring bars, PTFE, cylindrical, Diam. 6 mm, Length 20 mm	1 pz.
LLG09197525	Magnetic stirring bars, PTFE, cylindrical, Diam. 6 mm, Length 25 mm	
LLG09197530	Magnetic stirring bars, PTFE, cylindrical, Diam. 6 mm, Length 30 mm	1 pz.
LLG09197535	Magnetic stirring bars, PTFE, cylindrical, Diam. 6 mm, Length 35 mm	
LLG09197540	Magnetic stirring bars, PTFE, cylindrical, Diam. 8 mm, Length 40 mm	
LLG09197545	Magnetic stirring bars, PTFE, cylindrical, Diam. 8 mm, Length 45 mm	



Code	Description	Packaging
LLG09197550	Magnetic stirring bars, PTFE, cylindrical, Diam. 8 mm, Length 50 mm	
LLG09197560	Magnetic stirring bars, PTFE, cylindrical, Diam. 9 mm, Length 60 mm	
LLG09197570	Magnetic stirring bars, PTFE, cylindrical, Diam. 9 mm, Length 70 mm	
LLG09197580	Magnetic stirring bars, PTFE, cylindrical, Diam. 9 mm, Length 80 mm	



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

