

Technical Specification

PRODUCT : Carrez reagent zinc salt RS For agroalimentary analysis**CODE : 502700**

TEST	U.M.	SPECIFICATION
Density at 20°C	-	1.108 - 1.114
pH at 20°C	-	4.40 - 4.60
Zn(CH ₃ COOH)	: 219 g	
CH ₃ COOH	: 30 mL	
Water	: qsp 1000 mL	

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in