



Acétonitrile

Synonyms:

Méthyle cyanure, ACN,
CyanométhaneCH₃CN

Classification

Danger

Molecular weight

41.05 ONU

1648



CAS

75-05-8 Transport hazard
class:

3

EEC-N:

200-835-2 Packing

II

Acétonitrile

> RS - Pour HPLC - GOLD - Ultragradient

Grade Spécifique

Description	Clear liquid	-
Colour	<= 10	APHA
Identification	Positive	-
Miscible in acetone	Conform	-
Water miscibility	Conform	-
Diethyl ether miscibility	Conform	-
Miscibility in methanol	Conform	-



Acétonitrile

> RS - Pour HPLC - GOLD - Ultragradient

Grade Spécifique

Density at 20°C	0.781 - 0.785	-
Refractive index at 20°C	1.342 - 1.344	-
Boiling point	80.5 - 82.5	°C
Water (K.F.)	<= 100	ppm
Residue on evaporation	<= 2	ppm
Acidity	<= 0.0003	meq/g
Alcalinity	<= 0.0002	meq/g
Assay (GC)	>= 99.9	%
Fluorescence (quinine)	-	-
at 254 nm	<= 1	ppb
at 365 nm	<= 0.5	ppb
at 450 nm	<= 0.5	ppb
UV Absorbance	-	-
at 190 nm	<= 0.6	AU
at 200 nm	<= 0.03	AU
at 220 nm	<= 0.01	AU
at 254 nm	<= 0.005	AU
UV Transmittance	-	-
at 195 nm	>= 80	%
at 200 nm	>= 95	%



at 228 nm	>= 99	%
from 230 to 420 nm	>= 99	%
HPLC Gradient	-	-
at 210 nm	<= 1	mAU
at 254 nm	<= 0.2	mAU
Gradient baseline drift at 210 nm	<= 12	mAU
UV cut off	<= 190	nm

ID	Size	Packaging
000000000000412374		Flacon verre
000000000000412375		Bidon aluminium
000000000412371000		Flacon verre
412371000X6		Carton complet avec 6 x 1L flacon verre
000000000412372000		Flacon verre
412372000X4		Carton complet avec 4 x 2.5L flacon verre