



Dichlorométhane

Synonyms:
méthylène chlorure

CH₂Cl₂**Classification****Danger**

Molecular weight

84.93 ONU**1593**

CAS

75-09-2 Transport hazard
class:**6.1**

EEC-N:

200-838-9 Packing**III****Dichlorométhane****> RE - Puro - Stabilisé à l'éthanol****Grade Technique**

Description	Clear colourless liquid	-
Colour	<= 10	APHA
Density d ₂₀ /4	1.320 - 1.330	-
Refractive index at 20°C	1.4194 - 1.4294	-
Acidity (HCl)	<= 5	ppm
Residue on evaporation	<= 20	ppm
Water (K.F.)	<= 200	ppm



Assay (GC) (without stabilizer)	>= 99.8	%
Ethanol	~ 0.2	%

ID	Size	Packaging
000000000000337315	5 l	Fût métal
000000000000528370	25 l	Bidon plastique
000000000000528371	200 l	Fût métal
000000000000528372	5 l	Bidon plastique
000000000000528377	2.5 l	Flacon verre