



Ethyle acétate

Synonyms:

CH₃COOC₂H₅**Classification****Danger**

Molecular weight

88 ONU**1173**

CAS

141-78-6 Transport hazard
class:**3**

EEC-N:

205-500-4 Packing**II****Ethyle acétate**> **RE - Puro****Grade Technique**

Description	Clear, colourless liquid	-
Colour	<= 10	APHA
Identification (I.R.)	Positive	-
Density d _{20/20}	0.898 - 0.902	-
Refractive index at 20°C	1.3699 - 1.3739	-
Acidity (acetic acid)	<= 30	ppm
Residue on evaporation	<= 20	ppm



Water (K.F.)	<= 500	ppm
Assay (GC)	>= 99.8	%
Description	Clear, colourless liquid	-
Colour	<= 10	APHA
Identification (I.R.)	Positive	-
Density d20/4	0.898 - 0.902	-
Refractive index at 20°C	1.3699 - 1.3739	-
Acidity (acetic acid)	<= 30	ppm
Boiling point	76 - 78	°C
Water (K.F.)	<= 500	ppm
Assay (GC)	>= 99.8	%
Ethanol	<= 0.04	%

ID	Size	Packaging
000000000000508221	1 l	Flacon verre
000000000000508222	2.5 l	Flacon verre
000000000000528294	25 l	Fût métal
000000000000528295	5 l	Bidon plastique
000000000000528296	25 l	Bidon plastique
000000000000528297	200 l	Fût métal
000000000000528299	10 l	Fût métal