



## Glycérine 90% (28°Bé)

Synonyms:

1,2,3-propanetriol

CH<sub>2</sub>OHCHOHCH<sub>2</sub>OH

Classification

Danger

Molecular weight

92 ONU

CAS

56-81-5 Transport hazard class:

EEC-N:

200-289-5 Packing

Glycérine 90% (28°Bé)

> RE - Puro

Grade Technique

Description	Clear colourless liquid	-
Identification	Positive	-
Density at 20°C	1.220 - 1.233	-
Chlorinated compounds	<= 30	ppm
Heavy metals (as Pb)	<= 5	ppm
Residue on ignition	<= 100	ppm



**Glycérine 90% (28°Bé)**> **RE - Puro****Grade Technique**

Appearance of solution	Conform	-
Acidity or alcalinity	Conform	-
Refractive index at 20°C	1.4490 - 1.4550	-
Aldehydes	<= 10	ppm
Esters	Conform	-
Sugars	Conform	-
Chloride	<= 10	ppm
Water (K.F.)	11.5 - 16.5	%
Assay	83.5 - 88.5	%

ID	Size	Packaging
000000000000346131	1 l	Flacon verre
000000000000346132	2.5 l	Flacon verre
000000000000346134	35 kg	Fût plastique

**CARLO ERBA**  
REAGENTS

**CARLO ERBA Reagents S.A.S.**  
Chaussée du Vexin – 27106 Val de  
Reuil cedex  
N° TVA: FR 63391048824  
Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



*CARLO ERBA Reagents  
operates with a Certified  
Quality Management  
System*