



n-Hexane

Synonyms:

CH₃(CH₂)₄CH₃**Classification****Danger**

Molecular weight

86 ONU**1208**

CAS

110-54-3 Transport hazard
class:**3**

EEC-N:

203-777-6 Packing**II****n-Hexane****> RE - Puro****Grade Technique**

Description	Clear colourless liquid	-
Identification	Conform	-
Thiophene	Conform	-
Density d _{20/4}	0.655 - 0.665	-
Refractive index at 20°C	1.3699 - 1.3799	-
Boiling point	68.2 - 69.2	°C
Water (K.F.)	<= 100	ppm



n-Hexane > RE - Puro		Grade Technique
Residue on evaporation	<= 20	ppm
Acidity (Caproic acid)	<= 40	ppm
Assay (GC)	>= 95	%
Colour	<= 10	APHA
Aromatic compounds	<= 50	ppm
Total sulphur (S)	<= 5	ppm

ID	Size	Packaging
0000000000000339751	1 l	Flacon verre
0000000000000339752	5 l	Bidon plastique
0000000000000339755	2.5 l	Flacon verre
0000000000000339756	18 kg	Fût métal
0000000000000339758	25 l	Fût métal
0000000000000339759	200 l	Fût métal