

**n-Hexane 99%**

Synonyms:

CH₃(CH₂)₄CH₃**Classification****Danger**

Molecular weight

86.18 ONU**1208**

CAS

110-54-3 Transport hazard
class:**3**

EEC-N:

203-777-6 Packing**II****n-Hexane 99%****> RS - Pour GC-MS****Grade Specificque**

Description	Clear colourless liquid	-
Refractive index at 20°C	1.373 - 1.377	-
Water (K.F.)	<= 50	ppm
Residue on evaporation	<= 2	ppm
Colour	<= 10	APHA
Assay (GC)	>= 99.0	%
GC-MS. Individ. peak (hexadecane)	<= 2	°g/L



Ret. range n-undecane to n-tetracontane	-	-
(scanning area 30-600 amu)	-	-

ID	Size	Packaging
000000000000447212	1 l	Flacon verre