



Méthanol

Synonyms:

Alcool méthylique

CH3OH	Classification	Danger
Molecular weight	32 ONU	1230  
CAS	67-56-1 Transport hazard class:	3 
EEC-N:	200-659-6 Packing	II

Méthanol

> RS - RSE - Pour usage électronique

Grade Spécifique

Description	Clear liquid	-
Colour	<=10	APHA
Identification (I.R.)	Positive	-
Water miscibility	Conform	-
Readily carbonizable substances	Conform	-
Assay (GC)	>= 99.9	%
Resistivity	>= 0.5	MOhm.cm



Density at 20°C	0.791 - 0.793	-
Boiling point	64.1 - 65.1	°C
Water (K.F.)	<= 500	ppm
Residue on evaporation	<= 10	ppm
Refractive index at 20°C	1.328 - 1.330	-
Acidity (formic acid)	<= 15	ppm
Alcalinity (NH ₃)	<= 1	ppm
Ethanol	<= 200	ppm
Chloride	<= 0.2	ppm
Carbonyl compounds (as CO)	<= 5	ppm
Phosphate	<= 0.5	ppm
Heavy metals (as Pb)	<= 0.2	ppm
Subs. reducing KMnO ₄	<= 2.5	ppm
Total sulphur	<= 1	ppm
Ag	<= 0.2	ppm
As	<= 0.01	ppm
Al	<= 0.05	ppm
Au	<= 0.05	ppm
B	<= 0.01	ppm
Ba	<= 0.1	ppm



Be	<= 0.02	ppm
Bi	<= 0.02	ppm
Ca	<= 0.5	ppm
Cd	<= 0.005	ppm
Co	<= 0.01	ppm
Cr	<= 0.01	ppm
Cu	<= 0.01	ppm
Fe	<=0.1	ppm
Ga	<= 0.02	ppm
In	<= 0.02	ppm
K	<= 0.1	ppm
Li	<= 0.02	ppm
Mg	<= 0.1	ppm
Mn	<= 0.01	ppm
Mo	<= 0.05	ppm
Na	<= 0.5	ppm
Ni	<=0.01	ppm
Pb	<= 0.01	ppm
Pt	<= 0.05	ppm
Sb	<= 0.01	ppm



Sn	<= 0.02	ppm
Sr	<= 0.02	ppm
Ti	<= 0.05	ppm
Tl	<= 0.05	ppm
V	<= 0.05	ppm
Zr	<=0.05	ppm

ID	Size	Packaging
0000000000000414914	2.5 l	Flacon verre
0000000000000414917	1 l	Flacon verre