



Méthanol

Synonyms:

Alcool méthylique

CH ₃ OH	Classification	Danger
Molecular weight	32 ONU	1230  
CAS	67-56-1 Transport hazard class:	3 
EEC-N:	200-659-6 Packing	II

Méthanol

> RE - Puro

Grade Technique

Description	Clear colourless liquid	-
Identification (I.R.)	Positive	-
Density at 20°C	0.791 - 0.793	-
Refractive index at 20°C	1.3278 - 1.3298	-
Boiling point	64.1 - 65.1	°C
Residue on evaporation	<= 50	ppm
Acidity (formic acid)	<= 20	ppm



Water (K.F.)	<= 0.05	%
Subs. reducing KMnO4	<= 3	ppm
Assay (GC)	>= 99.9	%

ID	Size	Packaging
000000000000309000	200 l	Fût plastique
000000000000309001	2.5 l	Flacon verre
000000000000309002	25 l	Bidon plastique
000000000000309003	200 l	Fût métal
000000000000309004	1 l	Flacon verre
000000000000309008	10 l	Bidon plastique
000000000000309009	160 kg	Fût plastique
000000000000528101	5 l	Bidon plastique