



Xylène, mélange d'isomères

Synonyms:

C8H10	Classification	Danger
Molecular weight	106.17 ONU	1307  
CAS	Transport hazard class:	3 
EEC-N:	905-588-0 Packing	III

Xylène, mélange d'isomères> RPE - Pour analyse - ISO - ACS -
Reag.Ph.Eur. - Reag.USP**Grade Analytique**

Description	Clear liquid	-
Colour	<=10	APHA
Identification (I.R.)	Positive	-
Readily carbonizable substances	Conform	-
Density at 20°C	0.864 - 0.870	-
Refractive index at 20°C	1.4947 - 1.4987	-



**Xylène, mélange
d'isomères****> RPE - Pour analyse - ISO - ACS -
Reag.Ph.Eur. - Reag.USP****Grade Analytique**

Water (K.F.)	<= 200	ppm
Residue after evaporation	<= 20	ppm
Acidity (benzoic acid)	<= 14	ppm
Alcalinity (NH3)	<= 2	ppm
Benzene	<= 50	ppm
Ethylbenzene	<= 25	%
Thiophene	<= 1	ppm
Toluene	<= 1	%
Total sulphur	<= 3	ppm
Al	<= 0.5	ppm
Ba	<= 0.1	ppm
Ca	<= 0.5	ppm
Co	<= 0.02	ppm
Cr	<= 0.02	ppm
Cu	<= 0.02	ppm
Fe	<= 0.1	ppm
Mg	<= 0.1	ppm
Mn	<= 0.02	ppm
Ni	<= 0.02	ppm
Pb	<= 0.05	ppm



Xylène, mélange d'isomères		> RPE - Pour analyse - ISO - ACS - Reag.Ph.Eur. - Reag.USP	Grade Analytique
Sn	<= 0.1	ppm	
Zn	<= 0.05	ppm	
Assay (isomeric mixture)	>= 98.5	%	
Any single impurity	<= 0.5	%	

ID	Size	Packaging
000000000000492301	1 l	Flacon verre
492301X6	6 x 1 L	Carton complet avec 6 x 1L flacon verre
000000000000492303	23 kg	Fût métal
000000000000492304	170 kg	Fût métal
000000000000492305	5 l	Bidon plastique
000000000000492306	2.5 l	Flacon verre
492306X4	4 x 2,5 L	Carton complet avec 4 x 2.5L flacon verre