

Technical Specification

PRODUCT : PROPIONIC ACID RE - Pure**CODE : 306250****MOLECULAR WEIGHT :74.08****CAS N° : 79-09-4****FORMULA : C3H6O2**

TEST	U.M.	SPECIFICATION
Description	-	Yellow clear liquid
Identification	-	Positive
Density at 20°C	-	0.993 - 0.997
Refractive index at 20°C	-	1.3854 - 1.3894
Assay (acidimetric)	%	>= 98
Residue on evaporation	ppm	<= 100

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in