

Technical Specification

PRODUCT : 2-Phenylethanol RE Pure**CODE : 308730****MOLECULAR WEIGHT :****CAS N° :****FORMULA :**

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	Positive
Density at 20°C	-	1.017 - 1.023
Refractive index at 20°C	-	1.529 - 1.535
Assay (GC)	%	>= 98

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in