

# Technical Specification

**PRODUCT : Dichloromethane RE Pure stabilized with ethanol**

**CODE : 337310**

**MOLECULAR WEIGHT :84,93**

**CAS N° : 75-09-2**

**FORMULA : CH<sub>2</sub>Cl<sub>2</sub>**

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Colour	APHA	<= 10
Refractive index at 20°C	-	1.4194 - 1.4294
Density d <sub>20</sub> /4	-	1.320-1.330
Acidity (HCl)	ppm	<= 5
Residue on evaporation	ppm	<= 20
Assay (GC) (w ithout stabilizer)	%	>= 99.8
Water (K.F.)	ppm	<= 200
Ethanol	%	~ 0.2

Boiling point : 39.3 - 40.3 °C

completely miscible in Alcohol and Diethyl ether

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.

CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in