

Technical Specification

PRODUCT : Copper (II) chloride dihydrate RE - Pure

CODE : 364500

MOLECULAR WEIGHT :170,47

CAS N° : 10125-13-0

FORMULA : CuCl₂.2H₂O

TEST	U.M.	SPECIFICATION
Description	-	Blue crystals
Identification (I.R.)	-	Positive
pH sol. 5% at 20°C	-	3 - 3.8
Sulfate	ppm	<= 50
Nitrogenous compounds	%	<= 0.004
As	ppm	<= 1
Ca	ppm	<= 20
Fe	ppm	<= 10
K	ppm	<= 20
Mg	ppm	<= 20
Ni	ppm	<= 50
Pb	ppm	<= 40
Assay	%	>= 99

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in