

Technical Specification

PRODUCT : Paraffin white soft ERBApharm - According to pharmacopoeia : Ph.Eur.

CODE : 388400

CAS N° : 8009-03-8

TEST	U.M.	SPECIFICATION
Description	-	White mass
Identification	-	Positive
Drop point	°C	35 - 70
Substance aspect	Ph.Eur.	Conform*
Consistency	Ph.Eur.	60 - 300*
Acidity or alkalinity	Ph.Eur.	Conform*
Polycyclic aromatic hydrocarbons	ppm	<= 75*
Sulfated ashes	%	<= 0.05

CHARACTERS :

Appearance : white or almost white, translucent, soft unctuous mass, transparent in thin layers even after cooling to 0 °C.

Solubility : practically insoluble in water, slightly soluble in methylene chloride, practically insoluble in ethanol (96 per cent) and in glycerol.

This batch has been released by authorized personnel

* Parameters marked with * are tested periodically

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in