

Technical Specification

PRODUCT : TRIFLUOROACETIC ACID RS - For LC/MS

CODE : 411540

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	Positive
Water (K.F.)	%	<= 0.05
Chloride	ppm	<= 50
Fluoride	ppm	<= 50
Sulfate	ppm	<= 50
Assay (acidimetric)	%	>= 99.9
UV Absorbance	-	-
at 260 nm	AU	<= 0.9
at 270 nm	AU	<= 0.10
at 280 nm	AU	<= 0.08
at 300 nm	AU	<= 0.06
at 320 nm	AU	<= 0.05

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in