

# Technical Specification

**PRODUCT : Ethanol absolute anhydrous RS -For HPLC Isocratic grade**

**CODE : 412520**

**MOLECULAR WEIGHT :46,07**

**CAS N° : 64-17-5**

**FORMULA : C<sub>2</sub>H<sub>5</sub>OH**

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification (I.R.)	-	Positive
Density at 20°C	-	0.790 - 0.793
Refractive index at 20°C	-	1.3602 - 1.3622
Boiling point	°C	~ 78
Acidity or alcalinity	meq/g	<=0.0002
Water (K.F.)	%	<=0.05
Residue on evaporation	ppm	<=5
Assay (GC)	%	>=99.9
UV Transmittance	-	-
at 210 nm	%	>=25
at 220 nm	%	>=50
at 230 nm	%	>=75
at 240 nm	%	>= 80
at 250 nm	%	>=90
at 270 nm	%	>=94
at 290 nm	%	>=97

Filtered at 0.2 micron

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.  
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in