

Technical Specification

PRODUCT : HEXANE mixture of isomers RS - For HPLC - Isocratic Grade

CODE : 412630

MOLECULAR WEIGHT :86.18

CAS N° : 110-54-3

FORMULA : CH₃(CH₂)₄CH₃

TEST	U.M.	SPECIFICATION
Clear, colourless liq. appearance	-	Conform
Identification	-	Conform
Colour	APHA	<= 10
Refractive index at 20°C	-	1.375 - 1.383
Water (K.F.)	ppm	<= 150
Aromatic compounds	ppm	<= 5
Total sulphur (S)	ppm	<= 5
Residue on evaporation	ppm	<= 5
UV Transmittance	-	-
at 240 nm	%	>= 94
at 250 nm	%	>= 98
at 260 nm	%	>= 98

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in