

Technical Specification

PRODUCT : Chloroform RS PESTIPUR For pesticide analysis Stabilized with ethanol

CODE : 438650

MOLECULAR WEIGHT :119,38

CAS N° : 67-66-3

FORMULA : CHCl3

TEST	U.M.	SPECIFICATION
Description	-	Clear liquid
Colour	APHA	<= 10
Identification (I.R.)	-	Positive
Water (K.F.)	%	<= 0.01
Acidity (HCl)	ppm	<= 5
Residue on evaporation	ppm	<= 5
Assay (GC)	%	>= 99.9
GC-ECD.Individual peak (Lindane)	ng/L	<= 3
Ret.range 1,2,4 trichlorobenzene	-	-
to decachlorobiphenyle	-	-
GC-NPD.Individual peak (Ethylparathion)	ng/L	<= 3
Retention time Atrazin to Coumaphos	-	-

Stabilized with 0.4 to 1.0% Ethanol

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in