

Technical Specification

PRODUCT : Paraffin 56-58°C Erbaplast (without DMSO) RS For histology CE IVD**CODE : 467950**

TEST	U.M.	SPECIFICATION
Description	-	Solid w hite
Colour	-	min 25
Melting point	°C	56 - 60
Penetration at 25 °C	dmm	10 - 20
Viscosity at 100 °C	cst	3.5 - 5.0

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in