

Technical Specification

PRODUCT : Safranine T RS - For microscopy - C.I. 50420**CODE : 477230****MOLECULAR WEIGHT :350,85****CAS N° : 477-73-6****FORMULA : C20H19CIN4**

TEST	U.M.	SPECIFICATION
Description	-	Dark powder
Maximum lambda absorption (ETOH 50%)	nm	530 - 534

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in