

# Technical Specification

**PRODUCT : Tin (II) chloride dihydrate RPE - For analysis - ACS**

**CODE : 485000**

**MOLECULAR WEIGHT :225,63**

**CAS N° : 10025-69-1**

**FORMULA : SnCl<sub>2</sub>.2H<sub>2</sub>O**

TEST	U.M.	SPECIFICATION
Description	-	White crystals or almost white
Identification	-	Positive
Solubility in HCl	-	Conform
Sulfate	-	Conform
Ca	ppm	<= 50
Fe	ppm	<= 30
K	ppm	<= 50
Na	ppm	<= 100
Pb	ppm	<= 100
Assay (oxydimetric)	%	98.0 - 103.0

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.  
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in