

Technical Specification

PRODUCT : Carrez reagent potassium salt RS For agroalimentary analysis**CODE : 502710**

TEST	U.M.	SPECIFICATION
------	------	---------------

Density at 20°C	-	1.056 - 1.062
-----------------	---	---------------

K4Fe(CN)6 . 3H2O : 106 g

Water : qsp 1000 mL

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in