

Technical Specification

PRODUCT : Methanol RS For HPLC Isocratic Grade

CODE : 525100

MOLECULAR WEIGHT :32,04

CAS N° : 67-56-1

FORMULA : CH₃OH

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Assay (GC)	%	>= 99.9
Acidity (formic acid)	ppm	<=20
Density at 20°C	-	0.791 - 0.793
Refractive index at 20°C	-	1.3278 - 1.3298
Boiling point	°C	63.6 - 65.6
Residue on evaporation	ppm	<= 10
Water (K.F.)	ppm	<= 500
UV Transmittance	-	-
at 210 nm	%	>=20
at 235 nm	%	>= 80
at 260 nm	%	>= 98

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in