

Technical Specification

PRODUCT : Dichloromethane RE Pure stabilized with ethanol

CODE : 528370

MOLECULAR WEIGHT :84,93

CAS N° : 75-09-2

FORMULA : CH₂Cl₂

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Colour	APHA	<= 10
Density d20/4	-	1.320 - 1.330
Refractive index at 20°C	-	1.4194 - 1.4294
Acidity (HCl)	ppm	<= 5
Residue on evaporation	ppm	<= 20
Water (K.F.)	ppm	<= 200
Assay (GC) (w ithout stabilizer)	%	>= 99.8
Ethanol	%	~ 0.2

Boiling point : 39.3 - 40.3 °C

Completely miscible in Alcohol and Diethyl ether

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in