

Technical Specification

PRODUCT : Ethanol 95° denaturated RE - Pure - According to European regulation

CODE : 528770

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Colour	APHA	<= 10
Density at 20°C	-	0.804 - 0.812
Assay (alcohometric) at 20°C	% (V/V)	94.8 - 96.8

Regulation 2016/1867 of 20th October 2016

Denaturation by European process for 100L:

1 L of isopropyl alcohol

1 L of methylethylketone

1 g of Bitrex

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in