

Technical Specification

PRODUCT : HEXANE mixture of isomers RE - Pure

CODE : 528940

MOLECULAR WEIGHT :86,18

CAS N° : 110-54-3

FORMULA : CH₃(CH₂)₄CH₃

| TEST | U.M. | SPECIFICATION |
|--------------------------|------|-------------------------|
| Description | - | Clear colourless liquid |
| Identification | - | Positive |
| Density at 20°C | - | 0.665 - 0.675 |
| Refractive index at 20°C | - | 1.375 - 1.383 |
| Boiling point | °C | 63 - 70 |
| Water (K.F.) | ppm | <= 200 |
| Aromatic compounds | ppm | <= 30 |
| Residue on evaporation | ppm | <= 50 |
| Total sulphur (S) | ppm | <= 5 |

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in