

Technical Specification

PRODUCT : Dichloromethane RS ATRASOL For analysis of volatile traces Stabilized with ethanol

CODE : P02932E

MOLECULAR WEIGHT :84.93

CAS N° : 75-09-2

FORMULA : CH₂Cl₂

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless solution
Refractive index at 20°C	-	1.422 - 1.426
Water (K.F.)	ppm	<= 100
Residue on evaporation	ppm	<= 3
Colour	APHA	<= 10
Free acid (as HCl)	ppm	<= 5
Assay (GC) (w ithout stabilizer)	%	>= 99.95
Ethanol	% (w /w)	0.1 - 0.4
Carbon tetrachloride	ppm	<= 20
Chloroform	ppm	<= 40
GC-ECD.Individual peak (Lindane)	ng/L	<= 2
Ret.range 1,2,4 trichlorobenzene	-	-
to decachlorobiphenyle	-	-
GC-FID.Individual peak (n-hexadecane)	µg/L	<= 2
Ret. range n-undecane to n-tetracontane	-	-
PAH test (according to ISO 17993)	-	Passes test
Test against standard (each 16 PAH)	µg/L	<= 0.1

For traces analysis.

LIST OF STANDARD PAH :

- Naphthalene
- Acenaphtene
- Fluorene
- Phenanthrene
- Anthracene
- Fluoranthene
- Pyrene
- Benzo(a)anthracene
- Chrysene
- Benzo(b)fluorene
- Benzo(k)fluoranthene
- Benzo(a)pyrene
- Dibenzo(a,h)anthracene
- Benzo(g,h,i)perylene
- Indeno(1,2,3-c,d)pyrene
- Acenaphthylene

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in