

Technical Specification

PRODUCT : BORIC ACID INDICATOR

TEST	U.M.	SPECIFICATION
<p>Composition for 1L</p> <p>Boric Acid [CAS 10043-35-3] : 20g Methyl red solution (0.5 g/L in water): 10 mL Methylene blue solution (1.5 g/L in water): 2 mL Deionized water: up to 1L</p>	-	-

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application. CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of inappropriate use.

Date: 12/05/2026

Approved electronically, valid without a signature

**QUALITY LABORATORY
 TECHNICIAN
 GWENDOLINE HALOT
 VAL DE REUIL**