

Technical Specification

PRODUCT : SULFURIC ACID 69 % RS - For milk analysis

CODE : PS0893

TEST	U.M.	SPECIFICATION
H2SO4 content	%	68.0 - 70.0
Clear, colourless solution	-	Conform
Density d20/4	-	1.587- 1.611

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in