

# Technical Specification

**PRODUCT : SULFURIC ACID 62 % RS - For milk analysis**

**CODE : PS0894**

TEST	U.M.	SPECIFICATION
H2SO4 content	%	61.0 - 63.0
Clear, colourless solution	-	Conform
Density d20/4	-	1.509 - 1.531

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.  
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in