



## Miscela Etanolo 95° / Isopropanolo

Synonyms:

|                  | Classification          | Danger                                                                                                                                                                     |
|------------------|-------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Molecular weight | ONU                     | 1993   |
| CAS              | Transport hazard class: | 3                                                                                                                                                                          |
| EEC-N:           | Packing                 | II                                                                                                                                                                         |

### Miscela Etanolo 95° / Isopropanolo

&gt; RPE - Per analisi

Grado Analitico

|                                 |                         |         |
|---------------------------------|-------------------------|---------|
| Description                     | Clear colourless liquid | -       |
| Water (K.F.)                    | <= 5                    | %       |
| Density d20/4                   | 0.78 - 0.80             | -       |
| Ethanol content                 | >= 80                   | % (w/w) |
| Propan-2-ol content             | <= 15                   | % (w/w) |
| Assay (alcoholic) Total at 20°C | 94 - 96                 | % (V/V) |



Miscela Etanolo 95° /  
Isopropanolo

> RPE - Per analisi

Grado Analitico

ID

Size

Packaging

000000000000414551

5 l

Tanica di plastica



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



*CARLO ERBA Reagents  
operates with a Certified  
Quality Management  
System*