



Taq polimerasi

Solis Biodyne

FIREPoI DNA Polymerase:

Taq polimerasi robusta ed affidabile adatta a tutte le applicazioni di routine.

HOT FIREPoI DNA Polymerase:

Hot Start Taq Polimerasi con un'alta specificità e capace di ridurre gli eventi di primer dimers. Le performance elevate unite ad un alta riproducibilità dei differenti lotti di produzione la rendono adatta all'utilizzo anche nel campo diagnostico



Code	Description	Packaging
FC6010100500	FIREPol (regular DNA Polymerase) 500 units	500 units
FC6010101000	FIREPol (regular DNA Polymerase) 1000 units	1000 units
FC6010102000	FIREPol (regular DNA Polymerase) 2000 units	2000 units
FC6041100115	FIREPol 5x Master Mix 7.5 mM (250 rxn) 1 ml	1 ml
FC6041100125	FIREPol 5x Master Mix 12.5 mM (250 rxn) 1 ml	1 ml
FC6041200115	FIREPol 5x Master Mix 7.5 mM ready to load (250 rxn) 1 ml	1 ml
FC6041200125	FIREPol 5x Master Mix 12.5 mM ready to load (250 rxn) 1 ml	1 ml



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

