



## Taq polimerasi

## Solis Biodyne

**FIREPoI DNA Polymerase:**

Taq polimerasi robusta ed affidabile adatta a tutte le applicazioni di routine.

**HOT FIREPoI DNA Polymerase:**

Hot Start Taq Polimerasi con un'alta specificità e capace di ridurre gli eventi di primer dimers. Le performance elevate unite ad un alta riproducibilità dei differenti lotti di produzione la rendono adatta all'utilizzo anche nel campo diagnostico



| Code         | Description                                                | Packaging  |
|--------------|------------------------------------------------------------|------------|
| FC6010100500 | FIREPol (regular DNA Polymerase) 500 units                 | 500 units  |
| FC6010101000 | FIREPol (regular DNA Polymerase) 1000 units                | 1000 units |
| FC6010102000 | FIREPol (regular DNA Polymerase) 2000 units                | 2000 units |
| FC6041100115 | FIREPol 5x Master Mix 7.5 mM (250 rxn) 1 ml                | 1 ml       |
| FC6041100125 | FIREPol 5x Master Mix 12.5 mM (250 rxn) 1 ml               | 1 ml       |
| FC6041200115 | FIREPol 5x Master Mix 7.5 mM ready to load (250 rxn) 1 ml  | 1 ml       |
| FC6041200125 | FIREPol 5x Master Mix 12.5 mM ready to load (250 rxn) 1 ml | 1 ml       |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

