

**End-point PCR mix****Solis Biodyne**

Soluzioni flessibile disponibili sia in versione ready-to-load che con buffer a differenti concentrazioni di magnesio per fornire soluzione pronte all'uso ed adatte a tutte le esigenze sperimentali.

FIREPol Master Mix

Kit per l'amplificazione del DNA contenenti tutti i reagenti necessari ed adatte anche a templati difficili da amplificare in quanto ricchi inGC.

HOT FIREPol Blend Master mix

Kit contenente una miscela di HOT FIREPol (Hot Start) ed una polimerasi proof-reading per migliorare l'affidabilità della mix di amplificazione.



Code	Description	Packaging
FC6010200500	HOT FIREPol (HOT Start DNA Polymerase) 500 units	500 units
FC6010201000	HOT FIREPol (HOT Start DNA Polymerase) 1000 units	1000 units
FC6042500115	5x HOT FIREPol Blend Master Mix RTL 7.5 mM with BSA (250 rxn) 1 ml	1 ml
FC6042500120	5x HOT FIREPol Blend Master Mix RTL 10 mM with BSA (250 rxn) 1 ml	1 ml
FC6042500125	5x HOT FIREPol Blend Master Mix RTL 12.5 mM with BSA (250 rxn) 1 ml	1 ml
FC6042500130	5x HOT FIREPol Blend Master Mix RTL 15 mM with BSA (250 rxn) 1 ml	1 ml
FC6042502015	5x HOT FIREPol Blend Master Mix RTL 7.5 mM with BSA (5000rxn) 20 ml	20 ml
FC6042502020	5x HOT FIREPol Blend Master Mix RTL 10 mM with BSA (5000rxn) 20 ml	20 ml
FC6042502025	5x HOT FIREPol Blend Master Mix RTL 12.5 mM with BSA (5000rxn) 20 ml	20 ml
FC6042502030	5x HOT FIREPol Blend Master Mix RTL 15 mM with BSA (5000rxn) 20 ml	20 ml
FC6042510015	5x HOT FIREPol Blend Master Mix RTL 7.5 mM with BSA (25000rxn) 100 ml	100 ml
FC6042510020	5x HOT FIREPol Blend Master Mix RTL 10 mM with BSA (25000rxn) 100 ml	100 ml
FC6042510025	5x HOT FIREPol Blend Master Mix RTL 12.5 mM with BSA (25000rxn) 100 ml	100 ml
FC6042510030	5x HOT FIREPol Blend Master Mix RTL 15 mM with BSA (25000rxn) 100 ml	100 ml
FC6042700115	5x HOT FIREPol Blend Master Mix 7.5 mM with BSA (250rxn) 1 ml	1 ml
FC6042700120	5x HOT FIREPol Blend Master Mix 10 mM with BSA (250rxn) 1 ml	1 ml
FC6042700125	5x HOT FIREPol Blend Master Mix 12.5 mM with BSA (250rxn) 1 ml	1 ml
FC6042700130	5x HOT FIREPol Blend Master Mix 15 mM with BSA (250rxn) 1 ml	1 ml
FC6042702015	5x HOT FIREPol Blend Master Mix 7.5 mM with BSA (5000rxn) 20 ml	20 ml
FC6042702020	5x HOT FIREPol Blend Master Mix 10 mM with BSA (5000rxn) 20 ml	20 ml



Code	Description	Packaging
FC6042702025	5x HOT FIREPol Blend Master Mix 12.5 mM with BSA (5000rxn) 20 ml	20 ml
FC6042702030	5x HOT FIREPol Blend Master Mix 15 mM with BSA (5000rxn) 20 ml	20 ml
FC6042710015	5x HOT FIREPol Blend Master Mix 7.5 mM with BSA (25000rxn) 100 ml	100 ml
FC6042710020	5x HOT FIREPol Blend Master Mix 10 mM with BSA (25000rxn) 100 ml	100 ml
FC6042710025	5x HOT FIREPol Blend Master Mix 12.5 mM with BSA (25000rxn) 100 ml	100 ml
FC6042710030	5x HOT FIREPol Blend Master Mix 15 mM with BSA (25000rxn) 100 ml	100 ml
FC6043300115	5x HOT FIREPol GC Mmix	1 ml
FC60433001155	HOT FIREPol® GC Master Mix 5x1ml	
FC6043400120	5x HOT FIREPol« MultiPlex Mix with 10 mM MgCl ₂ (250rxn) 1 ml	1 ml
FC6043402020	5x HOT FIREPol« MultiPlex Mix with 10 mM MgCl ₂ (5000rxn) 20 ml	20 ml
FC6043600120	5x HOT FIREPol« MultiPlex Mix with 10 mM MgCl ₂ RTL (250rxn) 1 ml	1 ml
FC6043602020	5x HOT FIREPol« MultiPlex Mix with 10 mM MgCl ₂ RTL (5000rxn) 20 ml	20 ml
FC608010000110	HOT FIREPol® Multiplex qPCR Mix 10x1ml	
FC60801000015	HOT FIREPol® Multiplex qPCR Mix 5x1ml	





CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

