

**cDNA synthesis kit****Solis Biodyne****Kit per la retrotrascrizione disponibili nei formati****Synthesis MIX**

Miscela ready-to-load contenenti tutti i componenti necessari ad una retrotrascrizione.

Particolarmente adatti a tutte le applicazioni più critiche in cui la riduzione dei passaggi e step di dispensazione è indispensabile per limitare la variabilità degli esperimenti.

Synthesis KIT

Kit contenenti tutti i reagenti necessari in vial separate.

I synthesis Kit permettono la massima flessibilità e sono adatti a tutte le applicazioni di retrotrascrizione



Code	Description	Packaging
FC6061500050	FIREScript RT cDNA synth. Kit, 50 rxn 50 rxn	50 rxn
FC6061500200	FIREScript RT cDNA synth. Kit, 200 rxn 200 rxn	200 rxn
FC6061700100	FIREScript RT cDNA synth. Mix, NO primers, 100 rxn 100 rxn	100 rxn
FC6061700500	FIREScript RT cDNA synth. Mix, NO primers, 500 rxn 500 rxn	500 rxn
FC6061800100	FIREScript RT cDNA synth. Mix, oligo (dT), 100 rxn 100 rxn	100 rxn
FC6061800500	FIREScript RT cDNA synth. Mix, oligo (dT), 500 rxn 500 rxn	500 rxn
FC6061900100	FIREScript RT cDNA synth. Mix, random, 100 rxn 100 rxn	100 rxn
FC6061900500	FIREScript RT cDNA synth. Mix, random, 500 rxn 500 rxn	500 rxn
FC6062000100	FIREScript RT cDNA synth. Mix, oligo (dT)/random, 100 rxn 100 rxn	100 rxn
FC6062000500	FIREScript RT cDNA synth. Mix, oligo (dT)/random, 500 rxn 500 rxn	500 rxn
FC6063500050	SOLIScript RT cDNA synthesis KIT 50 rxn	50 rxn
FC6063500200	SOLIScript RT cDNA synthesis KIT 200 rxn	200 rxn
FC6063700100	SOLIScript RT cDNA synthesis MIX 100 rxn	100 rxn
FC6063700500	SOLIScript RT cDNA synthesis MIX 500 rxn	500 rxn





CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

