

**Maxima Probe/ROX qPCR Master Mix (2X)****Thermo Scientific**

qPCR master mix pronta all'uso - ottimizzata per SYBR Green usando protocolli di ciclizzazione standard - supplementata con ROX.



Code	Description	Packaging
FL50K0231	Maxima Probe/ROX qPCR Master Mix (2X)	200 reaz. (*)
FL50K0232	Maxima Probe/ROX qPCR Master Mix (2X)	1000 reac (*)
FL50K0233	Maxima Probe/ROX qPCR Master Mix (2X)	4000 reaz. (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

