

**T4 Polynucleotide Kinase (10 U/μL)****Thermo Scientific**

T4 polinucleotide chinasi (T4 PNK) - catalizza il trasferimento di gamma-fosfato da ATP a gruppi 5'-OH di oligonucleotidi, ss e dsDNA e RNA.



Code	Description	Packaging
FL5EK0031	T4 Polynucleotide Kinase (10 U/μL)	500 U (*)
FL5EK0032	T4 Polynucleotide Kinase (10 U/μL)	2500 U (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

