

**T4 DNA ligase****Thermo Scientific**

T4 DNA Ligasi - catalizza la formazione di legami fosfodiesterici tra 5'-P e 3'-OH di frammenti di dsDNA.



Code	Description	Packaging
FL5EL0011	T4 DNA ligase	1000 U (*)
FL5EL0012	T4 DNA ligase	5000 U (*)
FL5EL0013	T4 DNA ligase	HC, 5000 U (*)
FL5EL0014	T4 DNA ligase	200 U (*)
FL5EL0016	T4 DNA ligase	LC,1000 U (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

