



RiboLock RNase Inhibitor (40 U/ μ L)

Thermo Scientific

RiboLock RNase Inhibitor inibisce in modo specifico l'attività di RNasi A, B, e C, proteggendo l'RNA in diverse condizioni.



Code	Description	Packaging
FL5EO0381	RiboLock RNase Inhibitor (40 U/μL)	2500 U (*)
FL5EO0382	RiboLock RNase Inhibitor (40 U/μL)	10000 U (*)
FL5EO0384	RiboLock RNase Inhibitor (40 U/μL)	24x2500U



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

